

## HM 225.05

### Flow in a pipe bend



#### Learning objectives/experiments

- investigation of the pressure curve at a 90° pipe bend
- determination of the static pressure at 29 pressure measuring points
- representation of the pressure distribution

#### Description

- investigation of the pressure curve at a 90° pipe bend
- accessory for aerodynamics trainer HM 225

When laying pipes it is essential that they are adapted to the circumstances of their environment, which means the pipes will necessarily include deflections in the form of bends. Changing the direction of flow in a pipe changes the pressure conditions. The pressure curve during a change in the flow direction is investigated using the example of a 90° pipe bend.


The experimental unit HM 225.05, when used in the aerodynamics trainer HM 225, allows the measurement of the static pressure at 29 pressure measuring points along the pipe bend. The transparent pipe bend has a constant rectangular cross-section with 10 pressure measuring points each on the top and bottom. Pressure measuring points are located in the region of the curvature on both sides: four on the left side and five on the right side. The pressure measuring points are connected to the tube manometers in HM 225 via the hoses supplied. The static pressures can be read on the tube manometers.

To illustrate the pressure distribution, the static pressure at a measuring point is related to the maximum pressure. The graphical representation of the pressure curve shows a low pressure along the inner radius and an overpressure along the outer radius.


The experimental unit is attached to the HM 225 trainer, simply and precisely with quick release fasteners.

# HM 225.05

## Flow in a pipe bend



1 quick connector for connection to HM 225, 2 transparent pipe bend, 3 pressure measuring points, 4 tube manometers (HM 225)



Pressure distribution in the pipe bend:  $p^*$  related pressure change; orange: low pressure area, red: high pressure area

### Specification

- [1] determining the pressure conditions in flow through a pipe bend
- [2] measurement of static pressure at 29 pressure measuring points along the bend
- [3] 29 pressure measuring points: 4 on the left side and 5 on the right side, 10 on the top side and 10 on the bottom side
- [4] accessory for aerodynamics trainer HM 225
- [5] 16 tube manometers of HM 225 for displaying the pressures

### Technical data

90° pipe bend: cross-section 50x100mm  
29 pressure measuring points

LxWxH: 230x180x200mm  
Weight: approx. 4kg

### Scope of delivery

- 1 experimental unit
- 1 set of laboratory hoses
- 1 set of instructional material

# HM 225.05

## Flow in a pipe bend

Required accessories

070.22500      HM 225      Aerodynamics trainer