


Michelson interferometer - high resolution (Item No.: P2220900)



interference, wavelength, diffraction index, speed of light, phase, virtual light source


Overview

Short description

Interference, wavelength, diffraction index, speed of light, phase, virtual light source.

Principle

With the aid of two mirrors in a Michelson arrangement, light is brought to interference. While moving one of the mirrors, the alteration in the interference pattern is observed and the wavelength of the laser light determined.





Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

Equipment

Position No.	Material	Order No.	Quantity
1	He/Ne Laser, 5mW with holder	08701-00	1
2	Power supply for laser head 5 mW	08702-93	1
3	Interferometerplate w precision drive	08715-00	1
4	Optical base plate with rubberfeet	08700-00	1
5	Photoelement f. opt. base plt.	08734-00	1
6	Adjusting support 35 x 35 mm	08711-00	4
7	Surface mirror 30 x 30 mm	08711-01	4
8	Beam splitter 1/1,non polarizing	08741-00	1
9	Holder f. diaphr./beam splitter	08719-00	1
10	DMM, auto range, NiCr-Ni thermocouple	07123-00	1
11	Lensholder f. optical base plate	08723-00	1
12	Lens, mounted, f +20 mm	08018-01	1
13	Magnetic foot for optical base plate	08710-00	6
14	Screen, white, 150x150mm	09826-00	1
15	Measuring tape, $I = 2 m$	09936-00	1

Task

Construction of a Michelson interferometer using separate components. The interferometer is used to determine the wavelength of the laser light. The contrast function *K* is qualitatively recorded in order to determine the coherence length with it.

Set-up and procedure

In the following, the pairs of numbers in brackets refer to the co-ordinates on the optical base plate in accordance with Fig. 1. These co-ordinates are only a rough guideline. Perform the experimental set-up according to Fig. 1. The recommended set-up height (beam path height) is 130 mm.

- The lens L [1,7] must not be in position when making the initial adjustments. Before beginning with the adjustment, mount the extra plate P (fine adjustment drive on a plate) on the optical base plate according to Fig.1 ;whereby the co-ordinate lines should coincide as precisely as possible.
- When adjusting the beam path with the adjustable mirrors M₁ [1,8] and M₂ [1,4], the beam is aligned with the 4th y coordinate of the base plate.
- Adjust the mirror M3 [10,4], initially without the beam splitter BS [7,4], such that the reflected beam strikes the same point on mirror M₂ from which it previously originated.
- Now, place the beam splitter BS with its metallized side facing mirror M₂ in the beam path in such a manner that a partial beam strikes mirror M₃ unchanged and the other partial beam strikes mirror M₄ [7,1] perpendicularly along the 7th x coordinate of the base plate.
- The beam which is reflected by mirror **M**₄ must now be adjusted with the adjusting screws such that it strikes the same point on screen **SC** [7,6.5] as the partial beam that originated at mirror **M3** and was subsequently reflected by the beam splitter **BS** [7,4]. A slight flickering of the luminous points which have been made to coincide indicates nearly exact adjustment.
- By placing the lens L [1,7] in the beam path the luminous points are expanded.
- Now observe the interference patterns on screen SC (stripes, circles).
- By meticulously readjusting the mirrors M_3 and M_4 using the adjusting screws, one obtains concentric cireies.

On determining the wavelength of the laser light:

- To perform this measurement, the path distance between the mirror **M**₃ and the beam splitter **BS** must be changed. In the process, the position of mirror **M**₃ is altered using a lever arm (lever transmission ratio approx. 20:1) and a micrometer screw (2 turns correspond to 1 mm), and thus the optical path length of the light beam is also changed.
- On changing the optical path lengths, one sees changes in the centre of the interference rings from maxima to minima and visa versa. Whether the path length increases or decreases becomes apparent in the following: for decreasing path length,

Printed: 14/12/2017 15:54:10 | P2220900

the centre represents a source of maxima and minima; or for increasing path lengths it is a sink for the interference maxima and minima.

- According to the theory a change from minimum to minimum occurs when the optical path length $A \cdot d$ is changed by A, Le. in the set-up used the distance between the beam splitter **BS** and the mirror **M**₃ changes by A/2.
- To determine the wavelength of laser light, the changes in the distance between M₃ and BS are measured (by reading the initial and final values on the micrometer screw) and the number of changes from minimum to minimum (or maximum to maximum) are counted.

On recording the contrast function:

- In this case, the screen **SC** is replaced by a photo cell **PD** for the determination of the contrast function *K*. To ensure that the photocell does not measure the intensity across different maxima and minima of the circular interference fringes, reduce the size of the slotted diaphragm with black tape such that only a small aperture of approximately 1 mm² remains in the middle.
- For this part of the experiment, make the room as dark as possible to keep the dark current of the photocell as low as possible.
- To determine the contrast function, measure the intensities of minima and maxima at varying the optical separation of the mirrors. Change the separation using only mirror **M**₄. This mirror is only to be moved along the 7th x co-ordinate.
- Measure the distance between mirrors and beam splitter with a measuring tape. During the repositioning procedure, the mirror must be readjusted at each new position (if necessary, initially without the lens, see above) such that the interference fringes again become visible.
- To measure the intensities of minima and maxima, alter the position of the mirror M₃ slightly using the micrometer screw so that one can see which minimal and maximal voltage values can be measured with the multimeter (measuring range approx. 500 mV).
- The difference in optical path length between the two mirrors and the beam splitter should be varied between 0 and 10 cm: Le., when the distance from M₃ to the beam splitter is approximately 13 cm, mirror M₄ should be at position [7,2.25] at its minimum distance of approximately 8 cm from the beam splitter BS and at position [7,1.25] at its maximum separation of circa 13 cm from it.
- In the process, one must take into consideration that the larger the separation differences are, the smaller the radii of the circular interference fringes are. Consequently, at large separation differences the measurement of the maximum and minimum intensities is uncertain and, as a result of the relatively large diaphragm aperture, subject to large errors.

Theory and evaluation

If two waves having the same frequency *w* but different amplitudes and different phases are coincident at one location, they superimpose to:

 $E(t) = A_1 \cdot \sin(\omega t - arphi_1) + A_2 \cdot \sin(\omega t - arphi_2)$

The resulting wave can be described by the following:

 $E(t) = A_1 \cdot \sin(\omega t - arphi)$

with the amplitude

 $E^2(t) = A_1^2 + A_2^2 + 2 \cdot A_1 \cdot A_2 \cdot \cos \delta$ (1)

and the phase difference


 $\delta = \varphi_1 - \varphi_2$.

In a Michelson interferometer, light is split by a semitransparent glass plate into two partial beams (amplitude splitting), reflected by two mirrors, again brought to interference behind the glass plate (Fig. 2). Since only extensive luminous spots can exhibit circular interference fringes, the light beam is expandedlens **L**.




Robert-Bosch-Breite 10 D - 37079 Göttingen

Printed: 14/12/2017 15:54:10 | P2220900



If one replaces the real mirror M_3 between the laser and the glass plate by a virtual image M_3' , which is formed by reflection at the glass plate, a point **P** of the light source appears as the points **P**' and **P**'' of the virtuallight sources L_1 and L_2 .

Due to the different light paths, using the designations in Fig. 3 the phase difference is given by:



$$\delta = rac{2\pi}{\lambda} \cdot 2 \cdot d \cdot \cos heta$$
 (2)

 $\boldsymbol{\lambda}$ is the wavelength of the laser light used.

According to (1), the intensity distribution for $A_1 = A_2 = A$

$$I \!=\! E^2 \!=\! 4 \cdot A^2 \cos^2 rac{\delta}{2}$$
 (3)

Maxima thus occur when δ is equal to a multiple of 2 π , hence with (2)

$$2 \cdot d \cdot \cos \theta = m \cdot \lambda; m = 1, 2, \dots$$
 (4)

I. e., there exist circular fringes for selected fixed values of m, and d, since θ remains constant (see Fig. 3). If one alters the position of the movable mirror **M**₃ (cf. Fig.1) such that d, e.g., decreases, according to (4), the circular fringe diameter would also diminish since m is indeed defined for this ring. Thus, a ring disappears each time d is reduced by λ 2. At d = 0 the circular fringe pattern disappears. If the reflecting planes of mirrors **M**₃ and **M**₄ are not parallel in the sense of Fig. 3, one obtains curved fringes, which change into straight fringes at d = 0.

On determining the wave length:

To measure the wave length of the light, count the circular fringe changes while moving the mirror with the micrometer screw (transmission ratio approx. 20:1).

In the process, a shift of the mirror by 43.157 μ m is measured and N = 135(1) circular fringe changes are counted.

$$\lambda = rac{2 \cdot d}{N} = rac{2 \cdot 43.157}{135} \ \mu \mathrm{m}$$

PHYWE excellence in science

Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107 Printed: 14/12/2017 15:54:10 | P2220900



From these values, the wavelength of light $\lambda = 639(10)$ nm is obtained.

Temporal coherence and contrast function:

The temporal coherence - the coherence time and length - of a laser can be determined with the aid of a Michelson interferometer.

As a result of the different optical distances which the light traverses in the interferometer, the laser light which has been split into two beams undergoes a temporal retardation τ and is then caused to interfere with itself. The coherence time τ_c is the lag time at which the wave trains are still capable of interference, thus


 $au < au_c$. (5)

The coherence length is therefore:

 $l_c = c \cdot au_c$ (6)

where c : speed of light.

The laser's resonator determines the possible oscillation modes via its resonance condition $L_R = n \cdot \lambda/2 (n = 1, 2, 3...)$. In the process, only those frequencies which lie within the natural emission spectrum of the amplification medium and above the threshold of resonator losses occur (see Fig.4).



The contrast between the bright and dark circular fringes of the interference patterns is a measure of the interference capability of light. This can be determined by applying the autocorrelation function of light.

If $E(r,t) = A \cdot e^{i(kr - \omega t + \varphi)}$ is the complex electrical field vector of the light wave at location rand at time t, it follows that the intensity I (except for a constant factor) is:

$$I \,{=}\, (E \,{\cdot}\, E^*)$$
 (7)

where E^* is the conjugated complex vector of E and $\langle \rangle$ is a temporal average.

In our case the following results for the two waves E_1 and E_2 in the Michelson interferometer:

 $I_{
m res} = ((E_1+E_2)\cdot(E_1+E_2)^*) = I_1 + I_2 + 2\cdot Re(E_1\cdot E_2^*)$ (8)

In our experiment E_1 and E_2 are identical in the ideal case, except for the temporal shift τ , therefore:


 $E_2(t) = E_1(t+ au)$ (9)

$$\Gamma(au) = \left(E(t) \cdot E^{\,*}(t+ au)
ight)$$
 (10)

 $\Gamma(au)$ is the autocorrelation function or also the Self-coherence

function of light in this case.

For the resulting intensity, the following is thus obtained:



Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

Printed: 14/12/2017 15:54:10 | P2220900



 $I_{
m res} = 2I + 2Re(\Gamma(au))$ (11)

The contrast in the interference pattern is given by:

$$K\!=\!rac{I_{
m max}\,-\!I_{
m min}}{I_{
m max}\,+\!I_{
m min}}$$
 (12)

The standardised self-coherence function is the complex degree of self-coherence:

$$\gamma(au) = rac{\Gamma(au)}{\Gamma(0)}$$
 (13)

so that with (9) and (10) the following results:

$$E_1 = A_1 \cdot e^{i(k_1 x - \omega_1 au)}
onumber \ E_2 = A_2 \cdot e^{i(k_2 x - \omega_2 t)}$$
 (14)

and therefore

 $K = |\gamma(au)|$ (14)

For an ideal planar monochromatic wave in the x direction, the following contrast function would result:

$$K\!=\!|\gamma(au)|\!=\!1$$

for $E = A \cdot e^{i(kx - \omega t)}$ with $\gamma(\tau) = e^{-i\omega\tau}$ This means that the coherence time and length would be infinitely long in the ideal case for a single frequency. However, in reality, the coherence length is limited by the naturalline width (in gas lasers primarily additionally by Doppler broadening) of the spectral lines.

If a laser oscillates in two modes having the frequencies ω_1 and ω_2 with

 $egin{aligned} &I_{ ext{max}} = 2 \cdot |+2| \cdot |\gamma(au)| \ &I_{ ext{min}} = 2 \cdot |-2| \cdot |\gamma(au)| \end{aligned}$

the coherence function is the following:

$$egin{aligned} &\Gamma(au) = ((E_1(t) + E_2(t)) \cdot (E_1(t+ au) + E_2(t+ au))*) \ &= \Gamma_1(au) + \Gamma_2(au) \ &= |A_1|^2 e^{-i\omega_1 au} + |A_2|^{2e^{-i}\omega_2 au} \ \end{aligned}$$

and the degree of self-coherence is given by:

$$\gamma(au)=rac{|A_1|^2}{|A_1|^2+|A_2|^2}e^{-i\omega_1 au}+rac{|A_2|^2}{|A_1|^2+|A_2|^2}e^{i\omega_2 au}$$
 (16)

For simplification's sake A ₁and A ₂are equal. With $\omega_2 = \omega_1 + \Delta \omega$ we obtain the following results for the contrast of the interference grating:

$$\begin{split} K &= |\gamma(\tau)| \\ &= \sqrt{\frac{1}{2}(1 + \cos(\Delta \omega \tau))} \end{split}$$
(17)
$$&= |\cos(\frac{\Delta \omega \cdot \tau}{2})| \end{split}$$


If the 5-mW laser is used, the following frequency separation of the axial modes results for a resonator length L of approximately 30 cm:

 $\Delta\omega\,{=}\,2\pi\cdot\Delta f\,{=}\,2rac{\pi\cdot c}{2L}ig({=}\,3.1~{
m GHz}ig)$ (18)



Printed: 14/12/2017 15:54:10 | P2220900





The propagation delay time T results from the mirror shift d:

$$au = rac{2d}{c}$$
 (19)

With (18) and (19) the following results for the contrast function K:

$$K = \left|\cos\left(\frac{2\cdot\pi\cdot d}{2\cdot L}
ight)
ight| = \left|\cos\left(\frac{\pi\cdot d}{L}
ight)
ight|$$
 (20)

For the 5-mW laser the following is valid:

- according to (18), the mode separation Δf is obtained.
- with

$$K = rac{I_{ ext{max}} - I_{ ext{min}}}{I_{ ext{max}} + I_{ ext{min}}} = |\cos(rac{\Delta\omega\cdot au}{2})|$$


one obtains for the critical delay time au_c :

$$egin{aligned} K = 0 = |\cos(rac{\Delta \omega \cdot au_c}{2})| \,\, ext{or} \ rac{\pi}{2} = rac{2\pi \cdot \Delta f \cdot au_c}{2} o au_c pprox 1 \,\, ext{ns.} \end{aligned}$$

• This results in a minimum of the contrast function (20) at:

$$egin{aligned} o = |\cos(rac{\pi \cdot d}{L})| \,\, or \ rac{\pi}{2} = rac{\pi \cdot d}{L} &
ightarrow d = rac{L}{2} pprox 15.0 \,\, ext{cm}. \end{aligned}$$

The experimental data for the contrast function K as a function of the mirror shift d are given in Fig. 6.



It becomes apparent that the theoretical maximum is not reached; this can be due to several factors:

1. If the mirrors are inadequately parallel, the contrast function only reaches a smaller maximum value (thus, does not reach



Printed: 14/12/2017 15:54:10 | P2220900



1) and a larger minimum value (hence, does not drop to zero), see Fig. 7.

- 2. The division of the beam splitter is not ideal, Le. 50:50. Therefore, the derivation of the contrast function would have had to have been modified.
- 3. The 5-mW laser does not only oscillate in 2 modes, but rather the amplification is sufficient to allow the laser to oscillate in three axial modes: this shortens the coherence time. (A spectral analysis has established that a third mode is only possible to a considerably lesser degree than the two primary axial modes!)
- 4. The aperture in front of the photodiode was not made small enough. Consequently, it covers and averages different intensity regions.




Fig. 7: Theoretical contrast function K of 2-mode laser under ideal and real conditions.

No.	I_{\min} in mV	$I_{ m max}$ in mV	d in cm	K _{acc.} (12)
1	20	220	0.2	0.835
2	23	214.7	2.0	0.8065
3	67	200.5	5.3	0.4990

Table 1: Experimental data

